

Revolutionary way to accurately & easily micro-dose topical medications for humans & pets




-3-

Topi-CLICK.com

Standard Luer Lock Connection

REFILL Indicator

Innovative TRIAD Shape Allows Precision **Application** Control



1 CLICK = 0.05 mL 2 CLICKs = 0.10 mL 3 CLICKs = 0.15 mL



The Topi-CLICK Micro® Difference

Accurate & Consistent Dose™

- Topi-CLICK Micro® is the most accurate among those in a patient simulation test.
- Eliminates user-variability
- · Consistent automated dosing throughout treatment

Automated Metered Dosing

 Audible, Tactile and Visual CLICKs indicate measure of dose at each point of the unique triad-shaped base (One CLICK = 0.05 mL)

Targeted Topical Therapy

- Innovative triad-shaped dispenser allows for precision handling during application
- Applicator reduces the risk of hand transference to others
- · Offers applicator selections for customized dosing

Select one of 3 applicators



Universal Rounded Applicator¹



Angled **Applicator**



Pin Point **Applicator**™

Viscosity Assurance

- Accuracy NOT affected by varying formulations
- Accuracy NOT affected by varying bases

Topi-CLICK Micro® 9 mL

Dispense/CLICK: 0.05 mL

Dispenser Volume: Up to 9.75 mL Total Delivery Volume: Up to 9.45 mL

Approx. Residual: 0.15 mL

Approx. Priming: 0.15 or 3 CLICKs

Topi-CLICK Micro® 5 mL

Dispense/CLICK: 0.05 mL

Dispenser Volume: Up to 5.5 mL

Total Delivery Volume: Up to 5.35 mL

Approx. Residual: 0.15 mL

Approx. Priming: 0.15 or 3 CLICKs

Made of FDA grade materials | Phthalate free | blocks UV and is Recyclable 🖄 .

Topi-CLICK_MICRO_010122

www.Topi-CLICK.com

TO ORDER:

© 2022 TEAM Technologies

info.TopiCLICK@teamtech.com

Base

877-870-8448