




- ✓ Flexible Easy-press Buttons™ to eject full troches.
- ✓ CLICK locking mechanism eliminates accidental opening.
- ✓ Translucent tray allows for easy view of doses remaining.
- ✓ Available in 4 colors to aid in medication differentiation.

Flexible Easy-press Buttons™
make it easy to eject full troches.



Patient Use Instructions



STEP 1:

Remove troche tray from protective cover.



STEP 2:

Use Easy-press Buttons™ to eject troche lozenge. Dispense troche lozenge dosage as directed by your physician or pharmacist.



STEP 3:

Slide protective cover over tray to activate the CLICK locking mechanism.

Topi-CLICK® Troche-CLICK™ is Available in 4 Colors



Multiple Color Options make it easy for patients to distinguish between their medications.